



# REGENLOOP

ADVANCED GROWTH FACTOR CONCENTRATE (GFC) THERAPY KIT



Clinically Tested

100% Safe & Sterile

Dermatologist Recommended

Backed by Science



Next-Gen Technology

Scientifically Engineered

Powered by Growth Factors



Advanced Delivery System



Proven Results

Trusted by Experts

Designed for Excellence

High Efficacy, Zero Compromise

## Why you need this treatment...?



01

In today's stressful living, Alopecia (Hair Loss) is a common problem that affects men and women

02

A lot man factors contribute to it, including hormonal imbalances, stress, nutritional deficiency, genetic conditions etc.

03

Current scope of treatments include topical medications that have to be consumed indefinitely

04

Hair transplant surgeries too are not highly popular because of post surgery problems

05

Regenloop GFC Therapy offers you proven results in limited number of sessions as decided by your doctor



## What is Regenloop therapy




Regenloop GFC is pure, highly safe engineered from your own blood for superior results in hair loss


it helps with hair fall reduction, and make hair thicker and fuller in volume

Your blood has growth factors stored inside platelets (type of cell in our blood) which scientists have been able to extract as GFC in specially designed kit

# + How REGENLOOP GFC Therapy + benefits you?



## + What are the possible Side effects?





# INSTRUCTION FOR GFC PREPARATION

5ML

STEP 1



## COLLECTION OF BLOOD

Collect 4ml of blood in each tube


STEP 2



## ACTIVATION OF GFC

Mix by Inverting gently 8-10 Times and allow to Stand upright for 30 Minutes

STEP 3



## SEPARATION OF GFC

Centrifuge all the Tube at 3500 RPM for 10 Minutes

STEP 4



## COLLECTION OF GFC

Collection of GFC Serum from the tubes into the Syringes (up to 3ml per tube)

STEP 5



## APPLICATION

GFC is ready to Inject.

